

## USP54 Antibody

Catalog No: #46026

Package Size: #46026-1 50ul #46026-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	USP54 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	USP54 Antibody detects endogenous levels of total USP54
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human USP54
Target Name	USP54
Other Names	C10orf29
Accession No.	Swiss-Prot#: Q70EL1NCBI Gene ID:159195
Uniprot	Q70EL1
GeneID	159195;
Calculated MW	187kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

## Images



Western Blot analysis of HEK293 lysates, using primary antibody at 1:1000 dilution.

---

Note: This product is for in vitro research use only