

LRP1B Antibody

Catalog No: #46033

Package Size: #46033-1 50ul #46033-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	LRP1B Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	LRP1B Antibody detects endogenous levels of total LRP1B
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human LRP1B
Target Name	LRP1B
Other Names	LRP1B, LRPDIT, LRP-1B, LRP-deleted in tumors, LRP-DIT
Accession No.	Swiss-Prot#: Q9NZR2NCBI Gene ID:53353
Uniprot	Q9NZR2
GeneID	53353;
Calculated MW	515kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis of LRP1B expression in A375 (A), HeLa (B), H9C2 (C), NIH3T3 (D), K562 (E) whole cell lysates.

Note: This product is for in vitro research use only