

LRP4 Antibody

Catalog No: #46034

Package Size: #46034-1 50ul #46034-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	LRP4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	LRP4 Antibody detects endogenous levels of total LRP4
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human LRP4
Target Name	LRP4
Other Names	LRP4, CLSS, LRP-4, SOST2, KIAA0816, MEGF7
Accession No.	Swiss-Prot#: O75096NCBI Gene ID:4038
Uniprot	O75096
GeneID	4038;
Calculated MW	212kDa
Concentration	1.0mg/mL
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

Application Details

WB (1/500 - 1/2000)

Images



Western blot analysis of LRP4 expression in HepG2 (A), mouse lung (B) whole cell lysates.

Note: This product is for in vitro research use only