

MED31 Antibody

Catalog No: #46039

Package Size: #46039-1 50ul #46039-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MED31 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	MED31 Antibody detects endogenous levels of total MED31
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human MED31
Target Name	MED31
Other Names	MED31, CGI-125, HSOH1, Mediator complex subunit 31, Mediator complex subunit SOH1, Soh1
Accession No.	Swiss-Prot#: Q9Y3C7NCBI Gene ID:51003
Uniprot	Q9Y3C7
GeneID	51003;
Calculated MW	16 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Note: This product is for in vitro research use only