

## MID1IP1 Antibody

Catalog No: #46042

Package Size: #46042-1 50ul #46042-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	MID1IP1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	MID1IP1 Antibody detects endogenous levels of total MID1IP1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human MID1IP1
Target Name	MID1IP1
Other Names	MID1IP1, MIG12, S14R, Spot 14-R, Spot 14-related protein, Protein STRAIT11499, G12-like, MID1 interacting protein 1, Mid1-interacting protein 1, STRAIT11499, THRSPL
Accession No.	Swiss-Prot#: Q9NPA3NCBI Gene ID:58526
Uniprot	Q9NPA3
GeneID	58526;
Calculated MW	20 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

## Images



Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4o½ over night

---

Note: This product is for in vitro research use only