

MRPL23 Antibody

Catalog No: #46044

Package Size: #46044-1 50ul #46044-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	MRPL23 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	MRPL23 Antibody detects endogenous levels of total MRPL23
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human MRPL23
Target Name	MRPL23
Other Names	MRPL23, L23MRP, L23mt, MRP-L23, Ribosomal protein L23-like, RPL23L
Accession No.	Swiss-Prot#: Q16540NCBI Gene ID:6150
Uniprot	Q16540
GeneID	6150;
Calculated MW	18 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Note: This product is for in vitro research use only