

MOT7 Antibody

Catalog No: #46053

Package Size: #46053-1 50ul #46053-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	MOT7 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	MOT7 Antibody detects endogenous levels of total MOT7
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human MOT7
Target Name	MOT7
Other Names	SLC16A6, Breast cancer protein BCY5, MCT 6, MCT6, Monocarboxylate transporter 7, MCT 7, MCT7, Monocarboxylate transporter 6
Accession No.	Swiss-Prot#: O15403NCBI Gene ID:9120
Uniprot	O15403
GeneID	9120;
Calculated MW	57 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western Blot analysis of HEK293 lysates, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only