

MPHOSPH8 Antibody

Catalog No: #46055

Package Size: #46055-1 50ul #46055-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MPHOSPH8 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	MPHOSPH8 Antibody detects endogenous levels of total MPHOSPH8
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human MPHOSPH8
Target Name	MPHOSPH8
Other Names	MPP8
Accession No.	Swiss-Prot#: Q99549NCBI Gene ID:54737
Uniprot	Q99549
GeneID	54737;
Calculated MW	97 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western Blot analysis of HeLa lysates, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only