

FUNDC1 Polyclonal Antibody

Catalog No: #29805

Package Size: #29805-1 50ul #29805-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	FUNDC1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human Mouse
Immunogen Description	A synthetic peptide of human FUNDC1 (NP_776155.1).
Other Names	FUNDC1; FUN14 domain containing 1
Accession No.	Swiss-Prot#:Q8IVP5NCBI Gene ID:139341
Uniprot	Q8IVP5
GeneID	139341;
Calculated MW	17KDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using FUNDC1 antibody at 1:1000 dilution

Background

This gene encodes a protein with a FUN14 superfamily domain. The function of the encoded protein is not known.

Note: This product is for in vitro research use only