

HBEGF Polyclonal Antibody

Catalog No: #29821

Package Size: #29821-1 50ul #29821-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	HBEGF Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human HBEGF
Other Names	HBEGF; DTR; DTS; DTSF; HEGFL; heparin binding EGF like growth factor
Accession No.	Swiss-Prot#:Q99075NCBI Gene ID:1839
Uniprot	Q99075
GeneID	1839;
Calculated MW	23kDa
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.


Application Details

WB 1 μ O 500 - 1 μ O 2000IHC 1 μ O 50 - 1 μ O 200

Images



Western blot analysis of HBEGF expression in U87 (A), mouse skeletal muscle (B) whole cell lysates.



Immunohistochemical analysis of HBEGF staining in human colon cancer formalin fixed paraffin embedded tissue section. The section was pre-treated using heat mediated antigen retrieval with sodium citrate buffer (pH 6.0). The section was then incubated with the antibody at room temperature and detected using an HRP conjugated compact polymer system. DAB was used as the chromogen. The section was then counterstained with haematoxylin and mounted with DPX

Note: This product is for in vitro research use only