

## MYH3 Antibody

Catalog No: #46061



Package Size: #46061-1 50ul #46061-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	MYH3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	MYH3 Antibody detects endogenous levels of total MYH3
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human MYH3
Target Name	MYH3
Other Names	MYH3, HEMHC, SMHCE, MYHC-EMB, MYHSE1, Myosin heavy chain 3, Myosin-3
Accession No.	Swiss-Prot#: P11055NCBI Gene ID:4621
Uniprot	P11055
GeneID	4621;
Calculated MW	224kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only