

## MYO1A Antibody

Catalog No: #46063

Package Size: #46063-1 50ul #46063-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	MYO1A Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	MYO1A Antibody detects endogenous levels of total MYO1A
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human MYO1A
Target Name	MYO1A
Other Names	MYO1A, BBMI, BBM-I, DFNA48, MIHC, MYHL, Myosin I heavy chain, Brush border myosin I, Myosin IA, Unconventional myosin-Ia
Accession No.	Swiss-Prot#: Q9UBC5NCBI Gene ID:4640
Uniprot	Q9UBC5
GeneID	4640;
Calculated MW	118kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

## Images



Western Blot analysis of various cell lysis. Primary Antibody was diluted at 1:1000. Secondary antibody was diluted at 1:10000

---

Note: This product is for in vitro research use only