

SESN3 antibody

Catalog No: #23170



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


Description

Product Name	SESN3 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IF
Species Reactivity	Hu
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 21 and 266 of human SESN3
Target Name	SESN3
Accession No.	Swiss-Prot:P58005Gene ID:143686
Uniprot	P58005
GeneID	143686;
Concentration	1.36mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 57kd  
Western blotting: 1:500-1:3000  
Immunofluorescence: 1:100-1:200


Images



Sample (30 ug of whole cell lysate)  
A: Hela  
7.5% SDS PAGE  
Primary antibody diluted at 1: 1000



Merged with DNA probe



Immunofluorescence analysis of paraformaldehyde-fixed A549, using SESN3 antibody at 1: 200 dilution.

---

Note: This product is for in vitro research use only