

MT-ND4L Antibody

Catalog No: #46076

Package Size: #46076-1 50ul #46076-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	MT-ND4L Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	MT-ND4L Antibody detects endogenous levels of total MT-ND4L
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human MT-ND4L
Target Name	MT-ND4L
Other Names	MT-ND4L, MTND4L, NADH dehydrogenase subunit 4L, NADH dehydrogenase 4L, NADH4L, NAD4L, ND4L
Accession No.	Swiss-Prot#: P03901NCBI Gene ID:4539
Uniprot	P03901
GeneID	4539;
Calculated MW	11 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Note: This product is for in vitro research use only