

GPM6B Antibody

Catalog No: #46094

Package Size: #46094-1 50ul #46094-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GPM6B Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	GPM6B Antibody detects endogenous levels of total GPM6B
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human GPM6B
Target Name	GPM6B
Other Names	M6B
Accession No.	Swiss-Prot#: Q13491NCBI Gene ID:2824
Uniprot	Q13491
GeneID	2824;
Calculated MW	29 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Note: This product is for in vitro research use only