

P2RX1 Antibody

Catalog No: #46117



Package Size: #46117-1 50ul #46117-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	P2RX1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	P2RX1 Antibody detects endogenous levels of total P2RX1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human P2RX1
Target Name	P2RX1
Other Names	P2RX1, ATP receptor, p2X receptor 1, p2X1, p2X1 purinoceptor, p2X1 receptor, Purinergic receptor, p2RX1, p2X purinoceptor 1, p2X receptor, subunit 1
Accession No.	Swiss-Prot#: P51575NCBI Gene ID:5023
Uniprot	P51575
GeneID	5023;
Calculated MW	45 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only