

PEX11A Antibody

Catalog No: #46121

Package Size: #46121-1 50ul #46121-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PEX11A Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	PEX11A Antibody detects endogenous levels of total PEX11A
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human PEX11A
Target Name	PEX11A
Other Names	PEX11A, HsPEX11p, Peroxin-11A, PEX11, Protein PEX11 homolog alpha, PEX11-ALPHA, PMP28
Accession No.	Swiss-Prot#: O75192NCBI Gene ID:8800
Uniprot	O75192
GeneID	8800;
Calculated MW	28 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Note: This product is for in vitro research use only