

PAOX Antibody

Catalog No: #46122

Package Size: #46122-1 50ul #46122-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PAOX Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	PAOX Antibody detects endogenous levels of total PAOX
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human PAOX
Target Name	PAOX
Other Names	PAOX, PAO, Polyamine oxidase
Accession No.	Swiss-Prot#: Q6QHF9NCBI Gene ID:196743
Uniprot	Q6QHF9
GeneID	196743;
Calculated MW	70 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C


Application Details

Western blotting: 1:1000-3000


Images



Western Blot analysis of extracts from rat stomach, K562 cells, using PAOX Polyclonal Antibody. Secondary antibody was diluted at 1:20000



Immunohistochemical analysis of paraffin-embedded human-liver, antibody was diluted at 1:100



Immunohistochemical analysis of paraffin-embedded human-liver, antibody was diluted at 1:100

Note: This product is for in vitro research use only