

## PEX12 Antibody

Catalog No: #46123

Package Size: #46123-1 50ul #46123-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	PEX12 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	PEX12 Antibody detects endogenous levels of total PEX12
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human PEX12
Target Name	PEX12
Other Names	PEX12, PAF-3, PBD3A, Peroxisome assembly factor 3, Peroxin-12, Peroxisome assembly protein 12, PAF3, Peroxin 12
Accession No.	Swiss-Prot#: O00623NCBI Gene ID:5193
Uniprot	O00623
GeneID	5193;
Calculated MW	41 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

## Images



---

Note: This product is for in vitro research use only