

PIGC Antibody

Catalog No: #46126

Package Size: #46126-1 50ul #46126-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PIGC Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	PIGC Antibody detects endogenous levels of total PIGC
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human PIGC
Target Name	PIGC
Other Names	PIGC, GPI2, PIG-C
Accession No.	Swiss-Prot#: Q92535NCBI Gene ID:5279
Uniprot	Q92535
GeneID	5279;
Calculated MW	34 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis of lysates from Jarkat cells, primary antibody was diluted at 1:1000, 4o½ over night

Note: This product is for in vitro research use only