

PLSCR2 Antibody

Catalog No: #46134

Package Size: #46134-1 50ul #46134-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PLSCR2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	PLSCR2 Antibody detects endogenous levels of total PLSCR2
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human PLSCR2
Target Name	PLSCR2
Other Names	PLSCR2, PL scramblase 2, Phospholipid scramblase 2
Accession No.	Swiss-Prot#: Q9NRY7NCBI Gene ID:57047
Uniprot	Q9NRY7
GeneID	57047;
Calculated MW	34 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis of lysates from PC12 cells, primary antibody was diluted at 1:1000, 4o½ over night

Note: This product is for in vitro research use only