

PELO Polyclonal Antibody

Catalog No: #30607



Package Size: #30607-1 50ul #30607-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	PELO Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human PELO (NP_057030.3).
Other Names	PELO; CGI-17; PRO1770; protein pelota homolog
Accession No.	Swiss-Prot#:Q9BRX2NCBI Gene ID:53918
Uniprot	Q9BRX2
GeneID	53918;
Calculated MW	43kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:1000 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using PELO at 1:1000 dilution.

Background

This gene encodes a protein which contains a conserved nuclear localization signal. The encoded protein may have a role in spermatogenesis, cell cycle control, and in meiotic cell division.

Note: This product is for in vitro research use only