

KCNE1L Antibody

Catalog No: #46140

Package Size: #46140-1 50ul #46140-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	KCNE1L Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	KCNE1L Antibody detects endogenous levels of total KCNE1L
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human KCNE1L
Target Name	KCNE1L
Other Names	KCNE1L, AMMECR2, KCNE1-like, AMMECR2 protein, KCNE5
Accession No.	Swiss-Prot#: Q9UJ90NCBI Gene ID:23630
Uniprot	Q9UJ90
GeneID	23630;
Calculated MW	15 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4o½ over night

Note: This product is for in vitro research use only