

KCNH7 Antibody

Catalog No: #46142

Package Size: #46142-1 50ul #46142-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	KCNH7 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	KCNH7 Antibody detects endogenous levels of total KCNH7
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human KCNH7
Target Name	KCNH7
Other Names	KCNH7, Eag related protein 3, ERG-3, HERG3, ERG3, Kv11.3, HERG-3, Eag-related protein 3
Accession No.	Swiss-Prot#: Q9NS40NCBI Gene ID:90134
Uniprot	Q9NS40
GeneID	90134;
Calculated MW	135kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C


Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis of 1) Rat Brain Tissue, 2) Mouse Brain Tissue, 3) Hela with KV11.3 Rabbit pAb diluted at 1:2,000.



Immunohistochemical analysis of paraffin-embedded Mouse BrainTissue using Kv11.3 Rabbit pAb diluted at 1:200.

Note: This product is for in vitro research use only