

KCNS1 Antibody

Catalog No: #46143

Package Size: #46143-1 50ul #46143-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	KCNS1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	KCNS1 Antibody detects endogenous levels of total KCNS1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human KCNS1
Target Name	KCNS1
Other Names	KCNS1, HKv9.1, KV9.1
Accession No.	Swiss-Prot#: Q96KK3NCBI Gene ID:3787
Uniprot	Q96KK3
GeneID	3787;
Calculated MW	58 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western Blot analysis of HEK293 lysates, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only