

S100G Antibody

Catalog No: #46158

Package Size: #46158-1 50ul #46158-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	S100G Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	S100G Antibody detects endogenous levels of total S100G
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human S100G
Target Name	S100G
Other Names	S100G, CABP, S100 calcium-binding protein G, Protein S100-G, S100D, CABP9K, CALB3, Calbindin D9K, Calbindin-D9K, S100 calcium binding protein G
Accession No.	Swiss-Prot#: P29377NCBI Gene ID:795
Uniprot	P29377
GeneID	795;
Calculated MW	9 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis of S100G expression in rat spleen (A) whole cell lysates.

Note: This product is for in vitro research use only