

PCGF3 Polyclonal Antibody

Catalog No: #30931

Package Size: #30931-1 50ul #30931-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PCGF3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse
Immunogen Description	Recombinant fusion protein of human PCGF3 (NP_006306.2).
Other Names	PCGF3; DONG1; RNF3; RNF3A; polycomb group ring finger 3
Accession No.	Swiss-Prot#:Q3KNV8NCBI Gene ID:10336
Uniprot	Q3KNV8
GeneID	10336;
Calculated MW	28kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using PCGF3 at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit (RM00021). Exposure time: 30s.

Background

The protein encoded by this gene contains a C3HC4 type RING finger, which is a motif known to be involved in protein-protein interactions. The specific function of this protein has not yet been determined.

Note: This product is for in vitro research use only