

REPS2 Antibody

Catalog No: #46170



Package Size: #46170-1 50ul #46170-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	REPS2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	REPS2 Antibody detects endogenous levels of total REPS2
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human REPS2
Target Name	REPS2
Other Names	REPS2, Partner of RalBP1, POB1, RalBP1-interacting protein 2
Accession No.	Swiss-Prot#: Q8NFH8NCBI Gene ID:9185
Uniprot	Q8NFH8
GeneID	9185;
Calculated MW	72 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only