

RAB13 Antibody

Catalog No: #46176

Package Size: #46176-1 50ul #46176-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RAB13 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	RAB13 Antibody detects endogenous levels of total RAB13
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human RAB13
Target Name	RAB13
Other Names	RAB13, GIG4, RAS-associated protein RAB13, Ras-related protein Rab-13
Accession No.	Swiss-Prot#: P51153NCBI Gene ID:5872
Uniprot	P51153
GeneID	5872;
Calculated MW	23 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Note: This product is for in vitro research use only