

RAB15 Antibody

Catalog No: #46177

Package Size: #46177-1 50ul #46177-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	RAB15 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	RAB15 Antibody detects endogenous levels of total RAB15
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human RAB15
Target Name	RAB15
Other Names	RAB15, Ras-related protein Rab-15
Accession No.	Swiss-Prot#: P59190NCBI Gene ID:376267
Uniprot	P59190
GeneID	376267;
Calculated MW	24 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis of lysates from A431 cells, primary antibody was diluted at 1:1000, 4o½ over night

Note: This product is for in vitro research use only