

## RHG24 Antibody

Catalog No: #46203



Package Size: #46203-1 50ul #46203-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	RHG24 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	RHG24 Antibody detects endogenous levels of total RHG24
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human RHG24
Target Name	RHG24
Other Names	ARHGAP24, DKFZP564B1162, Filamin-A-associated RhoGAP, p73, p73RhoGAP, RCGAP72, RhoGAP of 73 kDa, Sarcoma antigen NY-SAR-88, FILGAP, RC-GAP72
Accession No.	Swiss-Prot#: Q8N264NCBI Gene ID:83478
Uniprot	Q8N264
GeneID	83478;
Calculated MW	84 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only