

## BCL2L15 antibody

Catalog No: #23184

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	BCL2L15 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC IF
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 1 and 130 of BCL2L15
Target Name	BCL2L15
Accession No.	Swiss-Prot:Q5TBC7Gene ID:440603
Uniprot	Q5TBC7
GeneID	440603;
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

## Application Details


Predicted MW: 18kd

Western blotting: 1:500-1:3000

Immunohistochemistry: 1:50-1:500

Immunofluorescence: 1:100-1:200

## Images




Sample (30 ug of whole cell lysate)


A: Raji

12% SDS PAGE


Primary antibody diluted at 1: 1000



Immunohistochemical analysis of paraffin-embedded A549 xenograft, using BCL2L15 antibody at 1: 500 dilution.



Merged with DNA probe



Immunofluorescence analysis of paraformaldehyde-fixed A549, using BCL2L15 antibody at 1: 200 dilution.

Note: This product is for in vitro research use only