

## MT-ND4 Polyclonal Antibody

Catalog No: #31776

Package Size: #31776-1 50ul #31776-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	MT-ND4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of mouse MT-ND4 (NP_904337.1).
Other Names	MT-ND4; MTND4; ND4; NADH dehydrogenase, subunit 4 (complex I)
Accession No.	Swiss-Prot#:_NCBI Gene ID:17719
GeneID	:17719
Calculated MW	51kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

## Application Details

WB □ 1:500 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using MT-ND4 at 1:1000 dilution.

## Background

Note: This product is for in vitro research use only