

SLC6A12 Antibody

Catalog No: #46245

Package Size: #46245-1 50ul #46245-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	SLC6A12 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Hu, MsoO Rt
Specificity	Recognizes endogenous levels of SLC6A12 protein.
Immunogen Type	Peptide
Immunogen Description	KLH-conjugated synthetic peptide region of human SLC6A12
Target Name	SLC6A12
Other Names	SLC6A12, BGT, BGT-1, BGT1, Gat2 mouse, Betaine transporter, GAT2
Accession No.	Swiss-Prot#: P48065NCBI Gene ID:6539
Uniprot	P48065
GeneID	6539;
Calculated MW	75 kD
Concentration	1.0mg/mL
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Application Details

WB 1:500 - 1:1000

Images



Western blot analysis of SLC6A12 expression in NIH3T3 (A), LO2(B), SGC7901 (C) whole cell lysates.

Note: This product is for in vitro research use only