

SLC6A13 Antibody

Catalog No: #46246



Package Size: #46246-1 50ul #46246-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SLC6A13 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	SLC6A13 Antibody detects endogenous levels of total SLC6A13
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human SLC6A13
Target Name	SLC6A13
Other Names	SLC6A13, GAT-2, GABA transport protein, GAT2
Accession No.	Swiss-Prot#: Q9NSD5NCBI Gene ID:6540
Uniprot	Q9NSD5
GeneID	6540;
Calculated MW	68 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only