

SLC6A20 Antibody

Catalog No: #46250

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SLC6A20 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	SLC6A20 Antibody detects endogenous levels of total SLC6A20
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human SLC6A20
Target Name	SLC6A20
Other Names	Solute carrier family 6 member 20B£B¬SIT1, XT3, XTRP3
Accession No.	Swiss-Prot#: Q9NP91NCBI Gene ID:54716
Uniprot	Q9NP91
GeneID	54716;
Calculated MW	66 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only