

SFRS16 Antibody

Catalog No: #46263

Package Size: #46263-1 50ul #46263-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	SFRS16 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	SFRS16 Antibody detects endogenous levels of total SFRS16
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human SFRS16
Target Name	SFRS16
Other Names	CLASRP, CLASP, SWAP2, SFRS16
Accession No.	Swiss-Prot#: Q8N2M8NCBI Gene ID:11129
Uniprot	Q8N2M8
GeneID	11129;
Calculated MW	77 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western Blot analysis of HEK293 lysates, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only