

SFRS18 Antibody

Catalog No: #46264

Package Size: #46264-1 50ul #46264-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	SFRS18 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	SFRS18 Antibody detects endogenous levels of total SFRS18
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human SFRS18
Target Name	SFRS18
Other Names	C6orf111, PNISR, SRRP130
Accession No.	Swiss-Prot#: Q8TF01NCBI Gene ID:25957
Uniprot	Q8TF01
GeneID	25957;
Calculated MW	93kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western Blot analysis of HeLa lysates, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only