

ZNF280B Antibody

Catalog No: #46266

Package Size: #46266-1 50ul #46266-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ZNF280B Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	ZNF280B Antibody detects endogenous levels of total ZNF280B
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human ZNF280B
Target Name	ZNF280B
Other Names	ZNF280B, 5OY11.1, 5'OY11.1, SUHW2, Zinc finger protein 279, ZNF632, Zinc finger protein 632, D87009.C22.3, Zinc finger protein 280B, ZNF279
Accession No.	Swiss-Prot#: Q86YH2NCBI Gene ID:140883
Uniprot	Q86YH2
GeneID	140883;
Calculated MW	62 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western Blot analysis of HEK293 lysates, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only