

SV2C Antibody

Catalog No: #46268

Package Size: #46268-1 50ul #46268-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SV2C Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	SV2C Antibody detects endogenous levels of total SV2C
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human SV2C
Target Name	SV2C
Other Names	KIAA1054
Accession No.	Swiss-Prot#: Q496J9NCBI Gene ID:22987
Uniprot	Q496J9
GeneID	22987;
Calculated MW	82 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C


Application Details

Western blotting: 1:1000-3000

Images



Western Blot analysis of extracts from Jurkat cells, using SV2C Polyclonal Antibody. Secondary antibody was diluted at 1:20000



Immunohistochemical analysis of paraffin-embedded rat-brain, antibody was diluted at 1:100

Note: This product is for in vitro research use only