

TIGD2 Antibody

Catalog No: #46273

Package Size: #46273-1 50ul #46273-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TIGD2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	TIGD2 Antibody detects endogenous levels of total TIGD2
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human TIGD2
Target Name	TIGD2
Other Names	Tigger transposable element-derived protein 2
Accession No.	Swiss-Prot#: Q4W5G0NCBI Gene ID:166815
Uniprot	Q4W5G0
GeneID	166815;
Calculated MW	60 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western Blot analysis of HEK293 lysates, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only