

ELOF1 Antibody

Catalog No: #46280

Package Size: #46280-1 50ul #46280-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ELOF1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	ELOF1 Antibody detects endogenous levels of total ELOF1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human ELOF1
Target Name	ELOF1
Other Names	ELOF1
Accession No.	Swiss-Prot#: P60002NCBI Gene ID:84337
Uniprot	P60002
GeneID	84337;
Calculated MW	9 kDa
Concentration	1.0mg/mL
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis of ELOF1 expression in HEK293T (A), SP2/O (B), H9C2 (C) whole cell lysates.

Note: This product is for in vitro research use only