

TCEAL2 Antibody

Catalog No: #46281

Package Size: #46281-1 50ul #46281-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TCEAL2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	TCEAL2 Antibody detects endogenous levels of total TCEAL2
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human TCEAL2
Target Name	TCEAL2
Other Names	TCEAL2, My048, TCEA-like protein 2, MY0876G05
Accession No.	Swiss-Prot#: Q9H3H9NCBI Gene ID:140597
Uniprot	Q9H3H9
GeneID	140597;
Calculated MW	26 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Note: This product is for in vitro research use only