

## TMC7 Antibody

Catalog No: #46286

Package Size: #46286-1 50ul #46286-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	TMC7 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	TMC7 Antibody detects endogenous levels of total TMC7
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human TMC7
Target Name	TMC7
Other Names	Transmembrane channel-like protein 7
Accession No.	Swiss-Prot#: Q7Z402NCBI Gene ID:79905
Uniprot	Q7Z402
GeneID	79905;
Calculated MW	84 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

## Images



Western blot analysis of TMC7 expression in THP1 (A), NIH3T3 (B), H9C2 (C) whole cell lysates.

---

Note: This product is for in vitro research use only