

USP24 Antibody

Catalog No: #46291

Package Size: #46291-1 50ul #46291-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	USP24 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	USP24 Antibody detects endogenous levels of total USP24
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human USP24
Target Name	USP24
Other Names	USP24, KIAA1057, Ubiquitin specific protease 24, Ubiquitin thioesterase 24, Ubiquitin thiolesterase 24, Deubiquitinating enzyme 24
Accession No.	Swiss-Prot#: Q9UPU5NCBI Gene ID:23358
Uniprot	Q9UPU5
GeneID	23358;
Calculated MW	294 kDa
Concentration	1.0mg/mL
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

Application Details

WB (1/500 - 1/2000)

Images



Western blot analysis of USP24 expression in HeLa (A), Jurkat (B), HT29 (C), mouse kidney (D) whole cell lysates.

Note: This product is for in vitro research use only