

TRIT1 Conjugated Antibody

Catalog No: #C27899



Package Size: #C27899-AF350 100ul #C27899-AF405 100ul #C27899-AF488 100ul
 #C27899-AF555 100ul #C27899-AF594 100ul #C27899-AF647 100ul
 #C27899-AF680 100ul #C27899-AF750 100ul #C27899-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	TRIT1 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	most applications
Species Reactivity	Hu
Immunogen Description	Recombinant fusion protein of human TRIT1 (NP_060116.2).
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	TRIT1; GRO1; IPPT; IPT; IPTase; MOD5; hGRO1; tRNA isopentenyltransferase 1
Accession No.	Swiss-Prot#:Q9H3H1NCBI Gene ID:54802
Uniprot	Q9H3H1
GeneID	54802;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	53kDa
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Background

This gene encodes a protein that is targeted to the mitochondrion and modifies transfer RNAs (tRNAs) by adding a dimethylallyl group onto the adenine at position 37. This modification is important for maintaining the correct reading frame during protein translation. This gene is considered a tumor suppressor and its expression can decrease cell growth. Alternative splicing results in multiple transcripts variants, most of which are likely non-functional.

Note: This product is for in vitro research use only