

USP34 Antibody

Catalog No: #46294

Package Size: #46294-1 50ul #46294-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	USP34 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	USP34 Antibody detects endogenous levels of total USP34
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human USP34
Target Name	USP34
Other Names	USP34, Deubiquitinating enzyme 34, KIAA0729, Ubiquitin specific protease 34, Ubiquitin thioesterase 34, Ubiquitin thiolesterase 34, KIAA0570
Accession No.	Swiss-Prot#: Q70CQ2NCBI Gene ID:9736
Uniprot	Q70CQ2
GeneID	9736;
Calculated MW	404 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western Blot analysis of lysate of mouse brain, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only