

USP49 Antibody

Catalog No: #46296

Package Size: #46296-1 50ul #46296-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	USP49 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	USP49 Antibody detects endogenous levels of total USP49
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human USP49
Target Name	USP49
Other Names	USP49, Deubiquitinating enzyme 49, Ubiquitin specific protease 49, Ubiquitin thiolesterase 49, Ubiquitin thioesterase 49
Accession No.	Swiss-Prot#: Q70CQ1NCBI Gene ID:25862
Uniprot	Q70CQ1
GeneID	25862;
Calculated MW	79 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western Blot analysis of HEK293 lysates, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only