

USP51 Antibody

Catalog No: #46297



Package Size: #46297-1 50ul #46297-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	USP51 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	USP51 Antibody detects endogenous levels of total USP51
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human USP51
Target Name	USP51
Other Names	Deubiquitinating enzyme 51 Ubiquitin thioesterase 51 Ubiquitin-specific-processing protease 51
Accession No.	Swiss-Prot#: Q70EK9NCBI Gene ID:158880
Uniprot	Q70EK9
GeneID	158880;
Calculated MW	80 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only