

## nrf2(Phospho-S40) Antibody

Catalog No: #12811

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	nrf2(Phospho-S40) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purification using immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	nrf2(Phospho-S40) Antibody detects endogenous levels of nrf2 only when phosphorylated at S40
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from human nrf2 around the phosphorylation site of S40.
Target Name	nrf2(Phospho-S40)
Modification	Phospho
Other Names	NFE2L2, NF-E2-related factor 2, NFE2-related factor 2, NRF2
Accession No.	Swiss-Prot#: Q16236NCBI Gene ID: 4780
Uniprot	Q16236
GeneID	4780;
Target Species	human
Calculated MW	90 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only